


Video Encounter Without Translation

Prior to initiating this workflow, please ensure that there is an active login with Doximity.

- Patient is checked in and roomed by MA
- Open the Rooming tab in patient chart
- Select **"Launch Doximity"** under Connect box



- Doximity window will open in browser with the phone number auto populated.


Caller ID: (510) 471-5907 • TVHC UC



Email a video invite >

- Select the video icon to initiate call
- Text message will be sent through the platform to the patient
- If using DAX Copilot, initiate visit and place phone near computer microphone for note generation
- A new window will automatically open for the visit
- It may take some time for the patient to connect through the platform

- If necessary, you can call the patient to facilitate connection to the video visit by clicking the box on the visit room screen.



- When visit is complete, end call